

Product Datasheet

Human DEFB103A ELISA Kit (orb405922)

Description	Human β -defensins 3 ELISA Kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Tested Applications	ELISA
Target	DEFB103A
Note	For research use only.
Sample Types	serum, plasma, urine, cell culture supernates, tissue homogenates
Uniprot ID	P81534
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, urine, cell culture supernates, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm