

Product Datasheet

Human UTS2 ELISA Kit orb405926

Description	Human Urotensin II ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	UTS2
Note	For research use only.
Sample Types	serum, saliva, cerebrospinal fluid (CSF), tissue homogenates
Uniprot ID	O95399
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, saliva, cerebrospinal fluid (CSF), tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm