

Product Datasheet

Human IFN gamma ELISA Kit orb405941

Description	Human Interferon gamma; ,IFN-gamma; ELISA Kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml
Tested Applications	ELISA
Target	IFNG
Note	For research use only.
Sample Types	Cell Culture Supernates
Uniprot ID	P01579
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, cell culture supernates, urine, tissue homogenates and cerebrospinal fluid (CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm