

Product Datasheet

Human GGT1 ELISA Kit orb405948

Description	Human gamma glutamyl transpeptidase,GGT ELISA Kit
Reactivity	Human
Range	0.78 mU/ml-50 mU/ml
Sensitivity	0.195mU/ml
Tested Applications	ELISA
Target	GGT1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm