

Product Datasheet

Human CA 27-29 ELISA Kit orb405953

Description	Human cancer antigen 27-29 (CA 27-29) ELISA Kit
Reactivity	Human
Range	3 U/ml-100 U/ml
Sensitivity	0.8 U/ml
Tested Applications	ELISA
Target	CA 27-29
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm