

Product Datasheet

Human IL11 ELISA Kit orb405960

Description	Human Interleukin 11,IL-11 ELISA KIT
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81pg/ml
Tested Applications	ELISA
Target	IL11
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, cell culture supernates, urine, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm