

Product Datasheet

Human IL15RA ELISA Kit orb405968

Description	Human Interleukin-15 receptor subunit alpha(IL15RA) ELISA kit
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.8 pg/ml
Tested Applications	ELISA
Target	IL15RA
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm