

Product Datasheet

Human IL16 ELISA Kit orb405969

Description	Human Interleukin 16,IL-16 ELISA KIT
Reactivity	Human
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml
Tested Applications	ELISA
Target	IL16
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm