

Product Datasheet

Human IRAK4 ELISA Kit orb405985

Description	Human interleukin-1 receptor-associated kinase 4(IRAK-4)ELISA Kit
Reactivity	Human
Range	6.25 pg/ml-400 pg/ml
Sensitivity	1.56 pg/ml
Tested Applications	ELISA
Target	IRAK4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm