

## Product Datasheet

### Human IL20 ELISA Kit (orb405989)

<b>Description</b>	Human Interleukin 20(IL-20)ELISA Kit...
<b>Reactivity</b>	Human
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Sensitivity</b>	0.078 ng/mL
<b>Target</b>	IL20
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Uniprot ID</b>	<a href="#">Q9NYY1</a>
<b>Expiration Date</b>	Please enquire.