

---

## Product Datasheet

### Human IL-23 ELISA Kit (orb405993)

<b>Description</b>	Human Interleukin 23,IL-23 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	6.25 pg/mL-400 pg/mL
<b>Sensitivity</b>	1.56 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL-23
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm