

## Product Datasheet

### Human IL24 ELISA Kit orb405996

<b>Description</b>	Human Interleukin 24(IL-24)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6 pg/mL-1000 pg/mL
<b>Sensitivity</b>	3.9 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL24
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Uniprot ID</b>	<b>Q13007</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm