

Product Datasheet

Human IL-28B ELISA Kit orb406001

Description	Human Interleukin 28B(IL-28B) ELISA kit
Reactivity	Human
Range	1.56 pg/mL-100 pg/mL
Sensitivity	0.39 pg/mL
Tested Applications	ELISA
Target	IL-28B
Note	For research use only.
Sample Types	serum, plasma, cell culture supernats, tissue homogenates
Uniprot ID	Q8IZI9
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernats, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm