

Product Datasheet

Human IL29 ELISA Kit (orb406002)

| | |
|----------------------------|--|
| Description | Human Interleukin 29(IL-29)ELISA Kit |
| Reactivity | Human |
| Range | 15.6 pg/mL-1000 pg/mL |
| Sensitivity | 3.9 pg/mL |
| Tested Applications | ELISA |
| Target | IL29 |
| Note | For research use only. |
| Sample Types | serum, plasma, cell culture supernates, urine, cerebrospinal fluid (CSF) |
| Uniprot ID | Q8IU54 |
| Expiration Date | Please enquire. |
| Application Notes | • Sample Type: serum, plasma, cell culture supernates, urine, cerebrospinal fluid (CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm |