
Product Datasheet

Human IL4 ELISA Kit (orb406014)

| | |
|----------------------------|---|
| Description | Human Interleukin 4,IL-4 ELISA KIT |
| Reactivity | Human |
| Range | 6.25 pg/mL-400 pg/mL |
| Sensitivity | 1.56 pg/mL |
| Tested Applications | ELISA |
| Target | IL4 |
| Note | For research use only. |
| Sample Types | serum, tissue homogenates, cell culture supernates, saliva |
| Uniprot ID | P05112 |
| Expiration Date | Please enquire. |
| Application Notes | • Sample Type: serum, tissue homogenates, cell culture supernates, saliva • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm |