

Product Datasheet

Human IL4 ELISA Kit orb406014

Description	Human Interleukin 4,IL-4 ELISA KIT
Reactivity	Human
Range	6.25 pg/ml -400 pg/ml
Sensitivity	1.56 pg/ml
Tested Applications	ELISA
Target	IL4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, tissue homogenates, cell culture supernates, saliva• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm