

Product Datasheet

Human IL5 ELISA Kit orb406016

Description	Human Interleukin 5,IL-5 ELISA KIT
Reactivity	Human
Range	0.78 pg/mL-50 pg/mL
Sensitivity	0.195 pg/mL
Tested Applications	ELISA
Target	IL5
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P05113
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm