

Product Datasheet

Human IL5 ELISA Kit orb406016

| | |
|----------------------------|--|
| Description | Human Interleukin 5,IL-5 ELISA KIT |
| Reactivity | Human |
| Range | 0.78 pg/ml-50 pg/ml |
| Sensitivity | 0.195 pg/ml |
| Tested Applications | ELISA |
| Target | IL5 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |