

Product Datasheet

Human IL6R ELISA Kit orb406018

Description	Human Interleukin 6 receptor(IL-6R)ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	IL6R
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P08887
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm