

Product Datasheet

Human LT-B4 ELISA Kit orb406022

Description	Human Leukotriene B4,LT-B4 ELISA Kit
Reactivity	Human
Range	1.25 ng/mL-80 ng/mL
Sensitivity	0.31 ng/mL
Tested Applications	ELISA
Target	LT-B4
Note	For research use only.
Sample Types	serum, urine, cell culture supernates, tissue homogenates
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, urine, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm