

Product Datasheet

Human bordetella pertussis antibody(IgM) ELISA Kit orb406033

Description	Human bordetella pertussis antibody(IgM)ELISA kit
Reactivity	Human
Range	Qualitative
Tested Applications	ELISA
Target	bordetella pertussis antibody(IgM)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm