

Product Datasheet

Human PKM2 ELISA Kit orb406055

Description	Human Pyruvate Kinase,M2-PK ELISA Kit
Reactivity	Human
Range	0.9 U/ml-60 U/ml
Sensitivity	0.225U/ml
Tested Applications	ELISA
Target	PKM2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm