

Product Datasheet

Human C1QB ELISA Kit orb406065

Description	Human Complement 1q,C1q ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	C1QB
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P02746
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm