

Product Datasheet

Human C1QC ELISA Kit [Out of Stock] orb406070

Description	Human Complement C1q subcomponent subunit C(C1QC) ELISA kit
Reactivity	Human
Range	9.38 ng/ml-600 ng/ml
Sensitivity	2.34ng/ml
Tested Applications	ELISA
Target	C1QC
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm