

Product Datasheet

Human C1QTNF2 ELISA Kit orb406073

Description	Human Complement C1q tumor necrosis factor-related protein 2(C1QTNF2) ELISA kit
Reactivity	Human
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml
Tested Applications	ELISA
Target	C1QTNF2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm