

Product Datasheet

Human Complement 4 ELISA Kit orb406078

Description	Human Complement 4,C4 ELISA Kit
Reactivity	Human
Range	0.41 ng/ml-300 ng/ml
Sensitivity	0.41 ng/ml
Tested Applications	ELISA
Target	Complement 4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm