
Product Datasheet

Human Complement 4 ELISA Kit (orb406078)

Description	Human Complement 4,C4 ELISA Kit
Reactivity	Human
Range	0.41 ng/mL-300 ng/mL
Sensitivity	0.41 ng/mL
Tested Applications	ELISA
Target	C4
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	POC0L4
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm