

Product Datasheet

Human Complement Component 9 ELISA Kit orb406080

Description	Human Complement Component 9, C9 ELISA Kit
Reactivity	Human
Range	10.97 pg/mL-8000 pg/mL
Sensitivity	10.97 pg/mL
Tested Applications	ELISA
Target	C9
Note	For research use only.
Sample Types	serum, plasma, urine
Uniprot ID	P02748
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, urine • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm