

Product Datasheet

Human C3d ELISA Kit orb406085

Description	Human Complement Fragment 3d,C3d ELISA Kit
Reactivity	Human
Range	0.156 µg/mL-10 µg/mL
Sensitivity	0.039 µg/mL
Tested Applications	ELISA
Target	C3d
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm