

Product Datasheet

Human C3a ELISA Kit orb406087

Description	Human Complement fragment 3a,C3a ELISA Kit
Reactivity	Human
Range	31.25 ng/mL-2000 ng/mL
Sensitivity	7.81 ng/mL
Tested Applications	ELISA
Target	C3a
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm