

Product Datasheet

Human C3b ELISA Kit orb406088

Description	Human Complement fragment 3b,C3b ELISA Kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Sensitivity	0.78 ng/mL
Tested Applications	ELISA
Target	C3b
Note	For research use only.
Sample Types	serum, urine, cerebrospinal fluid (CSF), saliva
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, urine, cerebrospinal fluid (CSF), saliva • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm