

Product Datasheet

Human LAMB1 ELISA Kit orb406102

Description	Human Laminin β 1(LN- β 1)ELISA Kit
Reactivity	Human
Range	9.38 ng/ml-600 ng/ml
Sensitivity	2.34 ng/ml.
Tested Applications	ELISA
Target	LAMB1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm