

Product Datasheet

Human ALP ELISA Kit orb406109

Description	Human Intestinal-type alkaline phosphatase,ALPI ELISA Kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Sensitivity	0.78 ng/mL
Tested Applications	ELISA
Target	ALPI
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09923
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm