

Product Datasheet

Human oxytocin ELISA Kit orb406145

Description	Human oxytocin ELISA Kit
Reactivity	Human
Range	8 µIU/ml-400 µIU/ml
Sensitivity	10 µIU/ml
Tested Applications	ELISA
Target	Oxytocin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm