

Product Datasheet

Human HSVI+II antibody (IgM) ELISA Kit orb406154

Description	Human herpes simplex virus I+II (HSVI+II) antibody (IgM) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	HSVI+II antibody (IgM)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm