

Product Datasheet

Human CETP ELISA Kit orb406160

Description	Human Cholesteryl ester transfer protein,CETP ELISA Kit
Reactivity	Human
Range	0.78 ng/ml-50 ng/ml
Sensitivity	0.195 ng/ml
Tested Applications	ELISA
Target	CETP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm