

Product Datasheet

Human Elastin ELISA Kit orb406165

Description	Human Elastin,ELN ELISA Kit
Reactivity	Human
Range	3.12 pg/ml-200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	Elastin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates, cell lysates, urine• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm