

## Product Datasheet

### Human PARK7 ELISA Kit orb406167

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Protein DJ-1 (PARK7) ELISA kit  |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 0.235 ng/mL-15 ng/mL  |
| <b>Sensitivity</b>         | 0.058 ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | PARK7   |
| <b>Note</b>                | For research use only.  |
| <b>Sample Types</b>        | serum, tissue homogenates, cell lysates   |
| <b>Uniprot ID</b>          | <a href="#">Q99497</a>  |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | • Sample Type: serum, tissue homogenates and cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm |