

Product Datasheet

Human FAM132A ELISA Kit orb406170

Description	Human Protein FAM132A(FAM132A) ELISA kit
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.8 pg/ml
Tested Applications	ELISA
Target	FAM132A
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm