

Product Datasheet

Human PRMT2 ELISA Kit orb406177

| | |
|----------------------------|---|
| Description | Human Protein arginine N-methyltransferase 2(PRMT2) ELISA kit |
| Reactivity | Human |
| Range | 37.5 pg/ml-2400 pg/ml |
| Sensitivity | 9.38 pg/ml. |
| Tested Applications | ELISA |
| Target | PRMT2 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |