

Product Datasheet

Human PRTN3 ELISA Kit orb406180

Description	Human proteinase 3(PR3) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Sensitivity	0.195 ng/mL
Tested Applications	ELISA
Target	PRTN3
Note	For research use only.
Sample Types	serum, tissue homogenates, plasma, urine
Uniprot ID	P24158
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, tissue homogenates, plasma and urine• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm