

Product Datasheet

Human BSA ELISA Kit orb406187

Description	Human bovine serum albumin (BSA) ELISA kit
Reactivity	Human
Range	0.312 µg/mL-80 µg/mL
Sensitivity	0.312 µg/mL
Tested Applications	ELISA
Target	BSA
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm