

Product Datasheet

Human LRP1 ELISA Kit orb406190

| | |
|----------------------------|--|
| Description | Human low-density lipoprotein-receptor-related protein 1,LRP-1 ELISA Kit |
| Reactivity | Human |
| Range | 0.625 ng/mL-40 ng/mL |
| Sensitivity | 0.156 ng/mL |
| Tested Applications | ELISA |
| Target | LRP1 |
| Note | For research use only. |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | Q07954 |
| Expiration Date | Please enquire. |
| Application Notes | • Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm |