

Product Datasheet

Human LRP4 ELISA Kit orb406191

Description	Human low-density lipoprotein-receptor-related protein 4(LRP-4)ELISA Kit
Reactivity	Human
Range	62.5 pg/ml-4000 pg/ml
Sensitivity	15.6 pg/ml
Tested Applications	ELISA
Target	LRP4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm