

## Product Datasheet

### Human RETNLB ELISA Kit (orb406198)

<b>Description</b>	Human Resistin-like Molecule $\beta$ (RELM- $\beta$ ) ELISA Kit...
<b>Reactivity</b>	Human
<b>Range</b>	47 pg/mL-3000 pg/mL
<b>Sensitivity</b>	11.7 pg/mL
<b>Target</b>	RETNLB
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Uniprot ID</b>	<a href="#">Q9BQ08</a>
<b>Expiration Date</b>	Please enquire.