

Product Datasheet

Human FAS ELISA Kit orb406203

Description	Human Factor-related Apoptosis,FAS ELISA Kit
Reactivity	Human
Range	1.25 ng/ml - 80 ng/ml
Sensitivity	0.31ng/ml
Tested Applications	ELISA
Target	FAS
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm