

Product Datasheet

Human AIFM1 ELISA Kit (orb406205)

Description	Human apoptosis inducing factor,AIF ELISA Kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Tested Applications	ELISA
Target	AIFM1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O95831
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm