

Product Datasheet

Human TERT ELISA Kit orb406211

Description	Human Telomerase reverse transcriptase(TERT) ELISA kit
Reactivity	Human
Range	6.25 ng/ml-400 ng/ml
Sensitivity	1.56 ng/ml
Tested Applications	ELISA
Target	TERT
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm