

## Product Datasheet

### Human TERT ELISA Kit orb406211

<b>Description</b>	Human Telomerase reverse transcriptase(TERT) ELISA kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	TERT
<b>Note</b>	For research use only.
<b>Uniprot ID</b>	<b>O14746</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm