

Product Datasheet

Human PTN ELISA Kit orb406223

Description	Human pleiotrophin,PTN ELISA Kit
Reactivity	Human
Range	0.24 ng/mL-15 ng/mL
Sensitivity	0.06 ng/mL
Tested Applications	ELISA
Target	PTN
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P21246
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm